AGN ENTERPRISES

MANUFACTURER & SUPPLIER
ALL KINDS OF LABORATORY EQUIPMENTS





paraffin tank.

It has independent heating zones that are controlled through a microprocessor. It has a large storage tank, the tank incorporates the filter which ensures continuos flow of paraffin from the tank. The flow of paraffin from the tank can be controlled manually or through foot switch. The heated operation platform and the forceps holder makes it easy to embed the tissues. The unit has a storage tank, preservation boxes, forceps holder, Paraffin nozzle whose Heating can be controlled through microprocessor and it can be programmed for the weekdays for both ON and OFF. Salient Features • Large paraffin tank suitable for high throughput lab ● Illuminated Work Platform ● Easy to read and operate control panel ● Optional Cold Spot ● Mold warmer and cassette bath ● It is programmable for weekdays and each heating zones temperature can be individually set which makes it easy for the user

Separate heating zones for storage tank, 2 preservation boxes, paraffin nozzle, forceps holder and work platform • Excess Paraffin collecting tray collects the excess paraffin from different heating zones

Read More

SKU:

Price: ₹300,000.00

Stock: instock

Categories: Histopathology

Product Description

It has independent heating zones that are controlled through a microprocessor. It has a large storage tank, the tank incorporates the filter which ensures continuos flow of paraffin from the

AGN ENTERPRISES



MANUFACTURER & SUPPLIER
ALL KINDS OF LABORATORY EQUIPMENTS



OFF. Salient Features ● Large parattin tank suitable for high throughput lab ● Illuminated work Platform ● Easy to read and operate control panel ● Optional Cold Spot ● Mold warmer and cassette bath ● It is programmable for weekdays and each heating zones temperature can be individually set which makes it easy for the user ● Separate heating zones for storage tank, 2 preservation boxes, paraffin nozzle, forceps holder and work platform ● Excess Paraffin collecting tray collects the excess paraffin from different heating zones