
# AGN ENTERPRISES







Polished Showcase)

**Read More** 

SKU:

**Price:** ₹950.00

**Stock:** instock

Categories: Geography & Surveying



#### **Product Description**

ROCK & MINERAL SET (PM= Paper Mounted , PS= Polished Showcase)

# AGN ENTERPR















(72x100CM)

**Read More** 

SKU:

**Price:** ₹3,800.00

Stock: instock

Categories: Geography & Surveying



## **Product Description**

(Any One Chart as Below) Educational Chart (72x100CM)

## **ROD & FRAME APPARATUS**

**Read More** 

SKU:

**Price:** ₹26,000.00

Stock: instock

**Categories:** Psychology



## AGN ENTERPR















SKU:

**Price:** ₹1,150.00

Stock: instock

**Categories:** Mathematics



#### **Product Description**

## **ROOT PRESSURE APPARATUS** (Borosilicate Glass)


Read More

SKU:

Price: ₹300.00

Stock: instock

**Categories:** Biology



## **Product Description**

# AGN ENTERPRISES

MANUFACTURER & SUPPLIER
ALL KINDS OF LABORATORY EQUIPMENTS



**Price:** ₹1,400.00

Stock: instock

**Categories:** Applied Mechanics



## **Product Description**

# ROTA ROD WITH DIGITAL TIMER Four Compartment



**Product Description** 

Read More

SKU:

**Price:** ₹16,000.00

Stock: instock

**Categories:** Pharmacy

## AGN ENTERPRI











Neau More

SKU:

**Price:** ₹6,000.00

Stock: instock

**Categories:** Pharmacy

#### **Product Description**





SKU:

**Price:** ₹14,000.00

Stock: instock

**Categories:** Pharmacy



#### **Product Description**

## **ROTA ROD WITH DIGITAL TIMER Two Compartment**



SKU:

**Price:** ₹13,000.00

Stock: instock

**Categories:** Pharmacy



### **Product Description**

## AGN ENTERP



ALL KINDS OF LABORATORY EQUIPMENTS













#### Read More

SKU:

**Price:** 

Stock: instock

Categories: Uncategorized

#### **Product Description**

**Specification :-** Type : Electronic lift Temperature range : Ambient to 220°C Heating power : 1500W. Operating voltage: 220 - 240 V at 50/60 HZ Display: LCD display Temperature deviation : at 2°C. Rotation speed : adjustment 20 to 280 RPM. **Features :-** - Cordless power supply allows emptying and refilling of water bath without unplugging. - LCD display indicates set and actual values of heating bath temperature, rotation speed, and lift height adjustment. -Auto lift from heating bath at the time of power failure - Accommodate flask size from 50 ml to 5000 ml. - Multifunction comb clip - A plastic coating on glass (P+G) - Condensate trap to trap external moisture formed on condenser - Provided with flask of 1000 ml each flask with flask joint SJ29/32 - Evaporation flask feeding via stopcock. - Rust free body material - Infrared bidirectional communication between heating bath and rotary evaporator Vacuum Pump: Operating voltage: 220 - 240 V at 50/60 HZ Membrane Pump: Two stage PTFE with window heads Flow rate: 1.8 m3/hr. Vacuum: 5 mbar **Vacuum Controller:** Display: Touch screen Operating mode: Manual, Timer, Continuous, Drying etc. Features: - - Rotary knob / button Allows in both ways - Built in solvent library (extended version: 60 Solvents). - Allows to add and

## AGN ENTERPR









continuously **Accessories:-** - Cora & plug - Warranty Cara - Instruction Manual



#### ROTARY MICROTOME

It is designed to enhance manual sectioning. It is based on cross roller guides with high precision screw which provides accuracy and stability for sectioning. This Microtome has been designed for effortless manual sectioning using counter balanced handwheel. The spacious section waste tray collects sectioning waste during trimming and sectioning. Salient

**Features** ● Based on cross roller guides ● Handwheel lockable at any position ● High precision specimen feed • Lateral adjustment of knife holder • Ergonomically located handwheel for coarse advance

Counter balanced handwheel Technical

**Specification** ● Section Thickness Setting Range: 0.5-60μm ● Setting Values: 0.5-2μm increment step 0.5μm ● 2-20μm increment step 1μm ● 20-60μm increment step 5μm • Minimum setting value :0.5μm Supplied complete as above with disposable blade holder assembly suitable for low & high profile blades

& sample pack of disposable blades Read More

SKU:

**Price:** ₹350,000.00

**Stock:** instock

Categories: Histopathology

#### **Product Description**

It is designed to enhance manual sectioning. It is based on cross roller guides with high precision screw which provides accuracy and stability for sectioning. This Microtome has been designed for effortless manual sectioning using counter balanced handwheel. The spacious

# AGN ENTERPRISES



MANUFACTURER & SUPPLIER
ALL KINDS OF LABORATORY EQUIPMENTS



Range: 0.5-60 $\mu$ m  $\bullet$  Setting values: 0.5- $2\mu$ m increment step  $0.5\mu$ m  $\bullet$  2- $20\mu$ m increment step  $1\mu$ m  $\bullet$  20-60 $\mu$ m increment step  $5\mu$ m  $\bullet$  Minimum setting value :0.5 $\mu$ m Supplied complete as above with disposable blade holder assembly suitable for low & high profile blades & sample pack of disposable blades